

Echopulse

Product Data 2013

General Overview

The Echopulse device is a high-precision device, specially designed and manufactured by Theraclion to treat benign tumors such as breast fibroadenomas and thyroid nodules, using High Intensity Focused Ultrasound (HIFU).

The main components are:

- Visualization and Treatment Unit (VTU)
- Integrated ultrasound scanner
- Arm and scanning motors
- Coupling and cooling system
- Patient motion detector
- Touch screen user interface

A patient bed and a set of cushions are necessary to perform the treatment. These elements are not provided by Theraclion.

A 4-Step Procedure

Preparation

- Device preparation and coupling and cooling set-up
- Patient installation
- Patient sedation/anaesthesia

Visualization

- Manual positioning of the treatment probe over the area to be treated and target visualization
- High-precision 3D electronic target centring

Planning

- Target outlining on the real-time ultrasound imaging
- Planning of HIFU target sites

Treatment

- Sequential treatment delivery (Pause & Pulse)
- Real-time view of HIFU pulses
- Continuous ultrasound monitoring



- The millimetric precision of HIFU delivering ensures efficacy and patients' safety.
- The liquid conductive cooling prevents skin heating.

General Features

Safety
<ul style="list-style-type: none"> ▪ No ionizing radiation ▪ Smart cooling and coupling technology ▪ Built-in diagnostic imaging ▪ Integrated and real-time ultrasound imaging ▪ Real-time patient monitoring and control (treatment can be stopped at all time) ▪ Patient motion detection ▪ Energy delivery control ▪ Built-in control safeties ▪ Emergency stop button ▪ Company certificated ISO 13485 and Echopulse CE marked for thyroid nodules and breast fibroadenomas
Efficacy
<ul style="list-style-type: none"> ▪ Alignment of the therapeutic beam with ultrasound imaging ▪ Pre-visualization of focused pulses ▪ Real-time view of HIFU pulses ▪ Millimetric treatment precision ▪ Coupling kit to optimize image quality
Cost-effective
<ul style="list-style-type: none"> ▪ Conscious sedation / local anaesthesia ▪ A 1h treatment procedure on average ▪ Minimal staff required




HIFU

The equipment is based on focused ultrasound technology. HIFU is a process that allows the delivery of a large amount of acoustic energy to a focal point. This acoustic energy is converted into heat (80°C-95°C) into the selected area to produce necrosis of targeted tissue.

Not approved for sale in the US.

Key Components

E-Pack



Before each treatment a new E-Pack must be installed on the device.

E-Pack replacement is secured through RFID management.

The cooling system consists of a cooler, two pumps and a disposable set of tubing and pouch containing 500ml of fluid.

It allows:

- Ultrasound waves to be transmitted without loss or distortion
- To protect the skin from the heat of HIFU waves

The pump circulates the fluid in a closed circuit between the pouch and the VTU. A temperature sensor located in the VTU ensures that the fluid in the VTU will be properly cooled. Internal pressure is adjusted accordingly.

Echopulse User Interface



The touch screen user interface is intuitive and easy to use. Mounted on an articulated arm, the screen is easily reachable and orientated by the operator.

Real-time imaging enables the operator to monitor the treatment and see its progress.

It is possible to suspend or modify the treatment parameters at any time.

The main functions of the user interface are:

- To guide the operator for kit installation
- To outline the key anatomical structures
- To determine parameters for imaging, positioning and energy
- To outline the target and select areas to be treated
- To launch and monitor the treatment in real-time

Visualization and Treatment Unit (VTU)



The VTU is composed of:




An ultrasound probe enabling real-time imaging of the region of interest.

To perform sagittal and transversal imaging of the target, the VTU rotates around its acoustical axis. The rotation of the VTU can be performed manually or controlled by the software during the planning phase.

An HIFU transducer is also embedded within the same treatment head, in-line with the imaging probe, in order to ensure optimal treatment safety and efficacy.

Schemes

- Ultrasound imaging
- Computers
- Coupling and cooling system



Technical Data & International standards

SYSTEM	PARAMETER		SPECIFICATION	
General Specifications	Reference		THC800154-D THC800154-E	
	Equipment designation		High Intensity Focused Ultrasound Therapy Equipment	
	Medical application		Ablation of tissues by external application	
	Classification 93/42/EEC		Equipment class IIb	
	Treatment duration		1 hour (average)	
	Total weight			
	Model THC800154-D		350 kg ± 20 kg	
	Model THC800154-E		370 kg ± 20 kg	
	Dimensions	Length	1150 mm	
		Width	600 mm ± 5 mm	
		Height	1300 to 1800 mm ± 20 mm	
	Power consumption		1 kW	
	Recommended preventive maintenance cycle		6 months	
Environmental conditions for transport and storage		4 °C to 65 °C		
Environmental conditions of operations		15 °C to 30 °C		
Arm	Type		6 degrees of freedom	
			Weight compensation	
			Electromagnetic friction brakes	
Electronic cabinet	Type		Mobile on 4 wheels	
Control computer	Monitor		19" Touch screen	
Epack	Reference		THC900800	
VTU	HIFU	Transducer	Piezoelectric	
		Type	Automatic depth compensation feature	
		Frequency	3MHz (central)	
		Typical maximal treatment depth	24 mm	
		Typical maximal acoustic power	125 W	
	Imaging probe	Type	128 elements linear array 38 mm wide	
		Default mode	B-mode	
		Frequency range	7.5 MHz - 12 MHz	
		Default image depth	5 cm	
		Typical frame rate	12 frames/s	
	Head			Motorized movements (R, x, y, z)
				xy range: ± 25 mm
				z range: ± 10 mm
		Robotic		
Main standards applied	IEC 60601 Ed.3.0 Medical Electrical Equipment			
	ISO14971 Medical Devices - Application Of Risk Management To Medical Devices			
	IEC 62304 Medical device software - Software life cycle processes			
	IEC 62366 Medical devices - Application of usability engineering to medical devices			
IP classification	IP 20			
Device Certification	CE marked for thyroid nodules and breast fibroadenomas			